



Q - I M P L A N T[®]

TRINON
TITANIUM

Back to the roots



Q-POWERPOINTS

- a gap free one piece titanium implant
- atraumatic one-phase concept
- immediate temporisation and immediate loading possibilities
- reduced treatment time for both patient and dentist
- easy to follow treatment sequence
- clearly arranged titanium instruments

THE Q-IMPLANT-CONCEPT

Q-Implant System was developed for easy and efficient use in dental practice. Designed to provide simplicity and clarity, we have created a system that removes many of the disadvantages of earlier implant systems. This single-phase enossal implant is made of titanium and there are no complex components.

Its self-cutting thread achieves a high primary stability which allows the immediate placement of a temporary crown. The transgingival healing makes a second operation unnecessary and the sand-blasted surface reduces the time taken for osseointegration. These features, together with the mirror-finished implant head, aid wound healing and reduce post-operative infections.

THE IMPLANT HEAD

The top of the implant was designed with a 7° cone. It can be individually grinded and can therefore be treated like a natural tooth. To avoid rotation, 4 symmetrical axial slots were placed on top.

They also guide the wrench you use to screw in the implant. This gives you the possibility to supply the implant with a temporal suprastructure. Before impression, the slots should be blocked out with composite. Using a dental dam, you may even grind the top during or right after surgery. Easy mucosa- and gum- management is possible due to the mirror-finish. Your patient leaves your office with an aesthetic and functional solution on the day of surgery.

THE THREAD

Q-Implants are constructed with a self-cutting thread. Gentle bone-management and atraumatic insertion is possible by using minimal force.

Even in the spongy bone of the upper jaw you reach a high primary stability. To reach this goal, we attached maximum importance in the development of a progressive designed thread with axial milling cuts, reducing stress and strain effects.

This spreads jaw pressure equally to the bone. The axial cuts also avoid rotation. Bone shavings can be harvested and re-used during osteogenesis.


THE SURFACE

For maximum osseointegration the surface was etched and blasted with aluminiumoxide.

Q-Implant surgery steps


Mucosal surgery:

Using a TRINON scalpel, the mucosa is opened by a crestal cut. A fold of skin is then built to uncover the osseous structures. Alternatively a gum punch can be used.




Initial indentation of the bone:

Use a rotating drill to create an indentation on the surface of the bone. Alternatively use a trocar drill to punch the mucosa if flapless surgery technique is applied.




Pilot drill:

Considering both of optimal prosthetic and enossal position of the implant, the right direction of the pilot drill is chosen.




Shaping drill:

The shaping drill is performed with the selected diameter and length.



Insertion:

Taken from the sterile cover the implant is gently inserted into the prepared bone cavern with the insertion wrench. It may be helpful to use the handwheel or the handwrench. By its self-cutting design the implant gains maximum primary stability.



Suture:

After reaching the implants final position, the mucosa is stitched up tightly. This is not necessary if a gum punch was used.




The suprastructure:

A temporal suprastructure is recommended during the healing phase by means of using a MultiCap+ or a Silicone Cap. The final prosthetic suprastructure follows.




Prosthetics Procedure


After accomplished woundhealing and without inflammation in mucosal area, the prosthetical process can start. In case the implants head are grinded the prosthetical forming is settled with conventional methods. Normally the slots in the head get blocked out with composite and the standard dental impression is executed. The laboratory model and the dental transformation are the same as providing a grinded tooth.




If the implants head are not grinded the use of MultiCap+ is recommended. This multifunctional cap is shifted onto the implants head and the impression made with an impression spoon.




The laboratory implants (3QLab / 4QLab / 5QLab) are now inserted into the Muticap+ which are sticking in the impression material for casting the model. Finally the technical laboratory work can be performed in usual manners.



The choice of a thermoplastic synthetic material gives the dental technician the possibility to use the MultiCap+ as a burn-out basis for a variety of dental applications. MultiCap+ has a uniform material strength due to its conical shape and constant thickness. Several laboratory steps are considerably shortened by using these prefabricated moulds. MultiCap+ burns out completely at 850°C (1562°F) without residues, which allows the use of a wide variety of alloys. Weight In gold approx. 0,8 g (12,35 gr.), which allows for no waste!



After the processing of the model all prosthetic possibilities are available. (e.g. ceramics full bridge)





The Q-Box



The Q-Box was designed for clarity and simplicity. It can be easily stacked into an autoclave and its sturdy construction ensures stability even following many sterilisation procedures. Every system component has its own unmistakable place making it easy for the operating assistant to locate the correct component quickly and effi-

ciently. The titanium bowl can be used for blending bone augmentation material as well as the temporary storage of operating-instruments.

| | Pilot-drill | Shaping-drill | Drill for cortical bone | Implant-holder | Insertion wrench | Mechanical insertion wrench | Hand wheel | Driver handle | Gum punch |
|---|-------------|----------------------|-------------------------|----------------|------------------|-----------------------------|-------------|---------------|----------------------------|
| Steribox complete with instruments | | | | | | | | | |
| QBOXC2 | PDQ1 | SDQ3 SDQ4 | SDHQ3 SDHQ4 | QGRIP | IWQL | IKQ3 | HWQ1 | DHQ1 | 3QPUNCH 4QPUNCH |
| <i>incl. rotating drill</i> | 2X | 1X | 1X | 1X | 1X | 1X | 1X | 1X | 1X |

Q-Tom

with the Q-Tom set the implantologist receives a bone-spreading and bone-condensing equipment completing the Q-Implant product range. These osteotomes made of titanium convince through their delicate design, maximum biocompatibility and elegant handling. The angulation of the convex working-ends support the use in frontal as well as in molar region. The Q-Tom set includes seven osteotomes beginning with a diameter of 1,8 mm.

The diameter rises by 0,2 mm per instrument to the maximum diameter of 3 mm. Through the minimal spreading of only 0,2 mm the operative effort and at the same time the risk of bone fractures in the spongious region is lowered. The excellent visibility of the lasermarking (length 8, 10, 12, and 14 mm) gives maximum safety for insertion depth.



| Complete set Osteotom Q-Tom <i>No. 1 – No. 7</i> | Osteotom Q-Tom No. 1 <i>(Bone Spreader)</i> | Osteotom Q-Tom No. 2 <i>(Bone Spreader)</i> | Osteotom Q-Tom No. 3 <i>(Bone Spreader)</i> | Osteotom Q-Tom No. 4 <i>(Bone Spreader)</i> | Osteotom Q-Tom No. 5 <i>(Bone Spreader)</i> | Osteotom Q-Tom No. 6 <i>(Bone Spreader)</i> | Osteotom Q-Tom No. 7 <i>(Bone Spreader)</i> |
|--|--|--|--|--|--|--|--|
| | ∅ 1,8 mm | ∅ 2,0 mm | ∅ 2,2 mm | ∅ 2,4 mm | ∅ 2,6 mm | ∅ 2,8 mm | ∅ 3,0 mm |
| QTOMC | QTOM1 | QTOM2 | QTOM3 | QTOM4 | QTOM5 | QTOM6 | QTOM7 |

In every case Q-Implant®

The Q-Implant product range suits all indications.

Such as single tooth treatment and bridges, increase of pillars in reduced set of teeth and in edentulous jaws.

It might be loaded immediately after tooth extraction and serves perfectly for immediate suprastructure, because of the high primary stability that is achieved by its very specific thread design.

Q-Implant, the implant system of the present and future !



Q - I M P L A N T ®

a popular and effective
one-phase concept
ideal for immediate temporisation
and immediate loading

Two different neck heights;
Standard (4 mm)
and Short (2 mm)

in diameter 2,5; 3,5; 3,9; 4,5 and
5,6 mm

and length 8, 10, 12 or 14 mm
available



G I P

an innovative one-phase concept
appropriate for strongly atrophied
but wide lower and upper jaw
ideal for immediate temporisation
and immediate loading

neck height 0,6 mm


in diameter 6,3 and 7,0 mm

and length 4, 5, 6, 7 mm available



Q³-IMPLANT®

one-phase implant with ballpoint-head
 excellent to be incorporated in existing overdentures
 or its implant supported new creation
 in diameter 3,5 mm and 4,5 mm
 and length 8, 10, 12 or 14 mm available



QK-IMPLANT®

two-phase implant with inner cone
 universal in use, fits all indications
 suitable for immediate and delayed loading and temporisation
 compatible to many other systems
 in diameter 4,0 and 5,6 mm
 and length 8, 10, 12 or 14 mm available




Q²-IMPLANT®

two-phase implant with external hexagon
 universal in use, fits all indications
 suitable for immediate and delayed loading and temporisation
 compatible to many other systems
 in diameter 3,5; 3,75 and 4,5 mm
 and length 8, 10, 12 or 14 mm available

Surgery

Implants



Q-Implant 2,5


| length | Standard | | | | | Short | | | |
|--------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|
| | Ø 2,5 mm Art. No. | Ø 3,5 mm Art. No. | Ø 3,9 mm Art. No. | Ø 4,5 mm Art. No. | Ø 5,6 mm Art. No. | Ø 3,5 mm Art. No. | Ø 3,9 mm Art. No. | Ø 4,5 mm Art. No. | Ø 5,6 mm Art. No. |
| 8 mm | 25Q08 | 3Q08 | 39Q08 | 4Q08 | 5Q08 | 3QS08 | 39QS08 | 4QS08 | 5QS08 |
| 10 mm | 25Q10 | 3Q10 | 39Q10 | 4Q10 | 5Q10 | 3QS10 | 39QS10 | 4QS10 | 5QS10 |
| 12 mm | 25Q12 | 3Q12 | 39Q12 | 4Q12 | 5Q12 | 3QS12 | 39QS12 | 4QS12 | 5QS12 |
| 14 mm | 25Q14 | 3Q14 | 39Q14 | 4Q14 | 5Q14 | 3QS14 | 39QS14 | 4QS14 | 5QS14 |



GIP-IMPLANT

GIP-Implant


| Length | Ø 6,3 mm Art. No. | Ø 7,0 mm Art. No. |
|--------|-------------------|-------------------|
| 4 mm | 63GIP104 | GIP4 |
| 5 mm | 63GIP105 | GIP5 |
| 6 mm | 63GIP106 | GIP6 |
| 7 mm | 63GIP107 | GIP7 |



Q³-IMPLANT®

Q³ Implant with matrix


| length | Ø 2,3 mm Art. No. | Ø 3,5 mm Art. No. | Ø 4,5 mm Art. No. |
|--------|-------------------|-------------------|-------------------|
| 8 mm | 23Q08 | 3Q308 | 45Q308 |
| 10 mm | 23Q10 | 3Q310 | 45Q310 |
| 12 mm | 23Q12 | 3Q312 | 45Q312 |
| 14 mm | 23Q14 | 3Q314 | 45Q314 |



QK-IMPLANT®

QK Implant with cover screw


| length | Ø 4,0 mm Art. No. | Ø 5,6 mm Art. No. |
|--------|-------------------|-------------------|
| 8 mm | 4QK08 | 56QK08 |
| 10 mm | 4QK10 | 56QK10 |
| 12 mm | 4QK12 | 56QK12 |
| 14 mm | 4QK14 | 56QK14 |



Q²-IMPLANT®













Q² Implant with cover screw

| length | Ø 3,5 mm Art. No. | Ø 3,75 mm Art. No. | Ø 4,5 mm Art. No. |
|--------|-------------------|--------------------|-------------------|
| 8 mm | 35Q208 | 4Q208 | 45Q208 |
| 10 mm | 35Q210 | 4Q210 | 45Q210 |
| 12 mm | 35Q212 | 4Q212 | 45Q212 |
| 14 mm | 35Q214 | 4Q214 | 45Q214 |



Instruments

| | | | |
|--|--|--|--|
| <p>Pilotdrill</p>  <p>PDQ1</p> | <p>Shaping drill</p>  <p>for implants</p> <p> \varnothing 2,5 mm \varnothing 3,5 mm \varnothing 3,9 mm SDQ25 SDQ3 SDQ39 \varnothing 4,5 mm \varnothing 5,6 mm SDQ4 SDQ5 </p> | <p>Drill for cortical bone</p>  <p>for implants</p> <p> \varnothing 3,5 mm \varnothing 3,9 mm \varnothing 4,5 mm \varnothing 5,6 mm SDHQ3 SDHQ39 SDHQ4 SDHQ5 </p> | <p>Insertion wrench</p>  <p>short IWQS long IWQL</p> |
|  <p>Q5108</p> | <p>Insertion hole drill</p>  <p>L= 4mm 5mm 6mm 7mm Q5104 Q5105 Q5106 Q5107</p> | <p>Insertion hole drill for cortical bone</p>  <p>L= 4mm 5mm 6mm 7mm Q5104K Q5105K Q5106K Q5107K</p> | |
|  <p>PDQ1</p> | <p>Shaping drill</p>  <p>for implants</p> <p> \varnothing 3,5 mm \varnothing 4,5 mm SDQ3 SDQ4 </p> | <p>Drill for cortical bone</p>  <p>for implants</p> <p> \varnothing 3,5 mm \varnothing 4,5 mm SDHQ3 SDHQ4 </p> | <p>Insertion wrench</p>  <p>short Q3123 long Q3124</p> |
| | <p>Shaping drill</p>  <p> \varnothing 4,0 mm \varnothing 5,6 mm Q4101 Q4102 </p> | <p>Drill for cortical bone</p>  <p>\varnothing 5,6 mm Q4103</p> |  <p>short Q2123 long Q2124</p> |
| | <p>Shaping drill</p>  <p> \varnothing 3,75 mm \varnothing 4,5 mm Q2101 Q2102 </p> | | |

| | | | |
|---|---|--|--|
| <p>Mechanical insertion wrench</p>  <p>IKQ3</p> | <p>Hand wheel</p>  <p>HWQ1</p> | <p>Driver handle</p>  <p>DHQ1</p> | <p>Gum punch</p>  <p> \varnothing 3 mm \varnothing 4 mm \varnothing 5 mm 3QPUNCH 4QPUNCH 5QPUNCH </p> |
| <p>Mechanical Insertion wrench</p>  <p>Q3120</p> | | | <p>Gum punch</p>  <p> \varnothing 3 mm \varnothing 4 mm 3QPUNCH 4QPUNCH </p> |
| <p>Mechanical insertion wrench</p>  <p><i>short</i> Q2120</p>  <p><i>medium</i> Q2121</p>  <p><i>long</i> Q2122</p> | | | <p>Screwdriver</p>  <p><i>short</i> 20 mm Q2113 (0,9 mm) Q2114 (1,25 mm)</p>  <p><i>long</i> 30 mm Q2113L (0,9 mm) Q2114L (1,25 mm)</p> <p>Handpiece for screwdriver</p>  <p>Q2112</p> |



Q-IMPLANT®





| | | | |
|--|---|--|--|
| <p>Laboratory analog</p> <p> \varnothing 2,5 mm \varnothing 3,5 mm \varnothing 3,9 mm \varnothing 4,5 mm \varnothing 5,6 mm 25QLAB 3QLAB 39QLAB 4QLAB 5QLAB </p> | <p>Transfer cap Multicap+*</p> <p>white caps for dentist impression</p> <p>QCAP QCAPS slot no Slot</p> <p>blue caps for dental technician only</p> <p>CAP5 CAP6 Slot no Slot</p> | <p>Silicone Cap</p> <p>fits all systems</p> <p>Q5001 - S Q5002 - M Q5003 - L</p> | <p>Easy Temporary Cap</p> <p>Acrylic cap</p> <p>\varnothing 4 mm - 10,5 mm</p> <p>Temporary Prosthesis made in 5 minutes. Please ask for detailed information.</p> |
| <p>Laboratory analogue for Q³-Implants**</p> <p>Q3210</p> | <p>Transfer cap for Q³-Implants**</p> <p>Q3310</p> | <p>Matrix for Q³-Implants**</p> <p>with O-Ring green, standard O-Ring green, standard O-Ring black, hard O-Ring blue, very soft</p> <p>Q3320 Q3321 Q3322 Q3323</p> | |
| <p>Gingivaformer</p> <p>2 mm 4 mm Q4003 Q4004</p> | <p>Cover screw</p> <p>small big Q4001 Q4002</p> | | |
| <p>Gingivaformer \varnothing 5 mm</p> <p>3 mm 4 mm 5 mm 7 mm Q2010 Q2011 Q2012 Q2013</p> | <p>Cover screw</p> <p>Q2001</p> | | |

* also for QK und Q2 with Q-Cone abutment















** also for QK und Q2 mit Q-Ball abutment

QK-Prosthetics







| | | | | |
|----------------------------|---|--|--|--|
| <i>Transfer/Laboratory</i> | Impression post | Laboratory analogue | Technician analogue | Fixing screw |
| | <i>with screw</i>  Q4010 (for 4QK) Q4029 (for 5,6 QK) |  Q4011 (for 4QK) Q4028 (for 5,6 QK) |  Q4041 |  Q4027 |

Abutments





| | | | |
|---|--|---|---|
| <i>Abutments screwed</i> | Abutment <i>screwed, straight</i>  Q4015 | Plastic cap <i>burn-out, screwed</i> <i>length 10 mm</i>  bridge Q4016 crown Q4017 | Occlusal screw  Q4018 |
| <i>Abutments cemented incl. fixing screw.</i> | Abutment <i>cemented</i>  <i>length 5,5 mm</i> Q4012 | Plastic cap <i>burn-out, cemented</i> <i>length 7 mm</i>  bridge Q4013 crown Q4014 | |
| <i>Abutments angled incl. fixing screw.</i> |  18° Q4022 |  18° Q4023 | |
| <i>Q-Cone Abutments</i> |  0 mm Q4019 |  2 mm Q4021 | <i>to be fixed with</i> IWQS, IWQL <i>Prosthetics like</i> Q-Implant Standard and Short |
| <i>Q-Ball Abutment</i> | Q-Ball Abutments with matrix | | |
| |  4 mm Q4026 |  3 mm Q4025 |  2 mm Q4024 |
| |  0 mm Q4020 |  | <i>to be fixed with</i> Q2110, Q2111 and Q2112 <i>Prosthetics like</i> Q ³ -Implant |

Q2-Prosthetics



| <i>Transfer/Laboratory</i> | Impression post <i>with screw</i>  Q2030 | Laboratory analogue  Q2040 | Technician analogue  Q2041 | Laboratory screw  Q2060 |
|----------------------------|--|---|---|--|
|----------------------------|--|---|---|--|

Abutments





| | | | | |
|---|--|--|--|--|
| <i>Straight and angled Abutments incl. fixing screw</i> |  | | | |
| <i>UCLA Abutment</i> |  <p>Fixing screw <i>for UCLA and all other abutments</i></p> <p><i>with screw</i> Q2058</p> <p>Q2056</p> | | | |
| <i>Q-Cone Abutments</i> |  <p><i>to be fixed with IWQS, IWQL</i></p> <p><i>Prosthetics like Q-Implant Standard and Short</i></p> | | | |
| <i>Q-Ball Abutment</i> | <p>Q-Ball Abutment with matrix</p>  <p>Q2063 Q2062 Q2061 Q2059</p> <p><i>to be fixed with Q2110, Q2111 and Q2112</i></p> <p><i>Prosthetics like Q³-Implant</i></p> | | | |

Para-Tubes – suitable for all implant systems




An intelligent solution for perfect parallelism. A safe insertion is guaranteed with the right distance and optimal axis.


Suitable for all implant systems.

| | | | | |
|--|---|---|--|--|
| Q-Para-Tubes-Set complete <i>(QPT08, QPT10, QPT12) with axis gauge and titanium box</i> QPTC | Depth and axis gauge  <i>with laser marking in 8, 10 und 12 mm</i> QBML | Drilling guide  <i>8 mm axis distance</i> QPT08 | Drilling guide  <i>10 mm axis distance</i> QPT10 | Drilling guide  <i>12 mm axis distance</i> QPT12 |
| | | | | |

Additional options

| | | | | | | | |
|---------------------------|---|---|---|---|---|--|--|
| Additional options | Trocar drill  QX26 | Rotating drill  QX27 | Drill extender  QX28 | Titanium Gaugeball <i>4 mm</i>  4QBALL | X-Ray Pattern  QXRAY | Titanium bowl  QX29 | Depth gauge  QPARAPIN |
| | Universal Torque Ratchet  <i>scale 1:2</i> QX30 | | | Guide Tubes  \varnothing 2,4 mm length 10 mm straight QX32  \varnothing 2,4 mm length 10 mm with shoulder QX34  \varnothing 2,4 mm length 6 mm straight QX33 | | | |
| | Scalpel Handle  <i>scale 1:3</i> QX31 | | |  \varnothing 3,4 mm length 6 mm straight QX35  \varnothing 3,4 mm length 6 mm with shoulder QX36 | | | |

Trinon product range for medicine




Q-BONE-GRAFTING-SET

- suitable for onlay-plastic and mesh supported augmentation
- bone screws in diameter 1,0 und 1,3 mm
- colour coding
- including titanium bowl for blending augmentation material or temporary storage of instruments




Q-MESH

- 3-dimensional formed titanium mesh
- applicable to totally atrophied maxilla
- reduced operation time
- easy to use
- individual modulation



Q-MULTITRACTOR-TYPE KARLSRUHE

- modular, vertical titanium distractor
- pre-implantological augmentation of mandibula and maxilla
- innovative pin-basis-plate
- high stability
- minimally invasive surgery
- simplified treatment, reduced operation time



BONE-PIN-SYSTEM

- titanium pins in length 3 mm and 5 mm
- for attachment of titanium mesh, -foils and membranes
- appropriate for 3-dimensional bone reconstruction
- titanium mesh 0,1 und 0,2 mm
- titanium foil 20 μ und 40 μ




OSTEOSYNTHESIS SYSTEM

- appropriate for maxillofacial-surgery
- titanium bonescrews and boneplates
- screw diameter 1,0 / 1,2 / 1,5 / 1,8 / 2,0 and 2,3 mm
- screwtop optional with inner crosshead or inner square
- large variety of plates with thickness 0,6 mm and 1,0 mm



SCALPELS

- scalpel blades
- disposable scalpels



MULTI-F

- external Ortho-Fixator made of titanium
- modular component system
- easy to use
- distraction, compression, dynamisation

HIGH QUALITY TITANIUM PRODUCTS SINCE 1993



M E D I C I N E



I N D U S T R Y



C O N S U M E R